

Breathe  
Easy

## PRO-MIST

### NEBULIZER COMPRESSOR

Optimized, Easy operation, Robust design

The Nidek Pro-Mist Nebuliser Compressor is a second generation, redesigned and heavy duty device that brings prompt relief to asthma patients including children. It is a home health care device which is used to administer medication in the form of a mist directly to the patient's airways. It is commonly used in treating asthma, cystic fibrosis, lung diseases, upper respiratory infection, and other respiratory diseases.

#### Product Highlight

- High nebulization rate ensures optimized treatment
- Easy operation and maintenance
- Small and robust design
- High capacity compressor



## Technical Specification

Electrical Requirements:	220-240Vac; 50Hz
Power Consumption:	160VA Max
Dimensions:	102mm (H) x 191mm (w) x 173mm (L), 4"(H) x 7.5" (W) x 6.81J (D)
Weight:	1.86 Kg. (4.1 Lb)
Sound Level:	56 dBA
Maximum Pressure:	2.41 bar (35 psi)
Operating Pressure:	0.76bar (11 psi)
Flow Range:	<9 lpm (0.318cfm)
Operation Humidity Temperature:	5°C (41°F) to 40°C (104°F) 30% to 95% RH
Storage Humidity Temperature:	- 40°C (-40°F) to 70°C (158°F) 10% to 1000% RH
Particle Size:	50~106 kpa (7.3~15.4 psi)

Note: Conditions may vary based on altitude barometric pressure and temperature. Specification subject to change without notice:



### **Nidek Medical India Private Limited**

**Head Office:** 238B A.J.C. Bose Road, Kolkata 700020,  
Tel: (033) 4017 7000 (30 Lines), Fax: 91 (033) 2289 4217

**Delhi:** 238 Ansal Chambers II, 6, Bhikaji Cama Palace  
New Delhi 110066, Tel: (011) 2616 5742, 93501 83222, 93135 77444

**Chennai:** AB-17, 5th Street, Anna Nagar, Chennai 600040  
Tel: (044) 4266 3662, 93808 06777, 93459 70044

**Mumbai:** 211 B, Damji-Shamji Industrial Estate, 9 LBS Marg, Kurla (W), Mumbai 400070  
Tel: (022) 4215 5441 / 93204 22266

**Hyderabad:** H.NO-3-41-66, Sri Laxmi Nilayam, South Block, West Marredpally  
Secunderabad 500026, Tel: 93468 32266

**Lucknow:** 4/11A, Ground Floor, Vibhav Khand, Gomti Nagar, Lucknow-226010  
Tel : 93361 20941

E-mail: [sales@nidekindia.com](mailto:sales@nidekindia.com), Website: [www.nidekindia.com](http://www.nidekindia.com)