

BREATHE EASY

Pulmo-Mist II Nebuliser

Twenty first century living has its built-in stress factors that add up to chronic health hazards. Lung diseases are on the rise with adults and children including infants getting affected. Today, millions of Indians are living with chronic lung diseases such as emphysema, asthma, and chronic bronchitis. Therefore, many asthma patients with common respiratory diseases are advised to use home medical equipment for treating these symptoms.

The Nidek Pulmo-Mist II Nebuliser is a second generation, redesigned and more safe and effective device that brings prompt relief to asthma patients including children.

It is a home health care device which is used to administer medication in the form of a mist directly to the patient's airways. It is commonly used in treating asthma, cystic fibrosis, lung diseases, upper respiratory infection, and other respiratory diseases.



PULMO-MIST II NEBULISER

Breathe Easy

Key Features

- It weighs only 1.4 Kgs - among the lightest in the industry
- Delivers medications quickly and effectively with a Free Flow rate of 10 lpm delivered by a powerful, reliable Piston Driven Compressor
- Styled better - compact design, rounded corners, integrated carrying handle, built-in nebuliser holder
- Built-in Safety - patient is protected by a built-in Compressor Relief valve which allows maximum 30 PSI delivery
- Durable - ABS injection moulded cabinet, damage-proof recessed air outlet, no breakable lids
- Feel-safe Performance - fully reliable patented motor design by the leader compressor manufacturer in the respiratory industry. Comes with a 1-year warranty.
- User-Friendly - easy to remove re-usable filter; fits any luggage or backpack; blends into any environment; single-button operation
- Economical - provides big savings in terms of power consumption and price

Specifications

Size	15.2cm x 10.1 cm x 17.8cm
Weight	1.4 kgs
Compressor Pressure	30 PSI (Relief Valve Setting)
Free Air Flow	10 lpm
Compressor Type	Piston
Capacity for medication	6 ml
Particle size	0.4 – 0.45 (miu)m
Nebulization rate	0.3 – 0.35 ml per minute
Electrical Requirements	230V, 50Hz - 0.7amps
Power Consumption	50-60 watts
Operation temperature	10 Celcius to 45 Celcius
Storage temperature	-25 Celcius to 70 Celcius
Humidity range	5 to 95% relative humidity
Warranty	1 Year



Nidek Medical India Private Limited

Head Office: 238B A.J.C. Bose Road, Kolkata 700 020
Tel: (033) 4017 7000 (30 Lines), Mob: 93395 66444

Delhi: 322, 3rd floor, Ansal Chamber II, 6, Bhikaji Cama Place, New Delhi 110 066
Tel: (011) 2616 5742, Mob: 93501 83222

Chennai: W-66, 2nd Street, Anna Nagar, Chennai 600 040
Tel: (044) 4266 3662, Mob: 93808 06777, 93459 70044

Mumbai: 211 B, Damji-Shamji Industrial Estate, 9 LBS Marg, Kurla (W), Mumbai 400 070
Tel: (022) 4215 5441 / 80809 14400

Hyderabad: H.NO-3-41-66, Sri Laxmi Nilayam, South Block, West Marredpally, Secunderabad 500 026
Tel: (040) 4015 0243, 93468 32266

Lucknow: 1/44 Vijayant Khand Gomti Nagar, Lucknow 226 010
Tel: 93361 20941

Chandigarh: 2nd Floor, SCO# 355, Sector 32-D, Chandigarh 160 031
Tel: 93562 78899

E-mail: sales@nidekindia.com, Website: www.nidekindia.com

