

Cardiac Monitoring

NIDEK ECG Machine ECG-701

Portable, Precise, User-friendly

Nidek India's Electrocardiograph machines are one of the basic diagnostic tools for cardiologists and general physicians.

They can be used to perform routine diagnostic checks. The machines are simple to use and precise, thus, making the diagnosis easy to perform and analyze.

Nidek's ECG machines are available in several models suited to the physicians and patient's needs and requirements.

Product descriptions appear in the following pages.

Key Features

- Single channel ECG machine with interpretation
- Colour display
- Manual / Auto / RR / Store Modes (optional)
- Built in rechargeable battery
- Adapt to 110-230V, 50 / 60Hz AC power supply
- Built in RS232 interface (ECGNET software optional)



Specifications

Display	LCD 3.5"
Paper	50mm X 20m roll / fold paper
Dimension	285 X 265 X 72 Unit (mm)

Standard Specifications

Lead	Standard 12 leads
Acquisition	Simultaneous 12 leads
Resolution	12bit / 1000Hz
Operation mode	Manual / Auto / RR / Store Modes
Filter	AC, EMG Filter
Drift filter	Anti-Drift System: 0.15 Hz (Adaptive)
CMRR	$\geq 120\text{dB}$, with AC Filter
Input circuit	Floating, Protection against defibrillator effect
Input circuit current	$\leq 0.05 \mu\text{A}$
Input Impedance	$> 50 \text{M}\Omega$
Patient current leakage	$< 10 \mu\text{A}$
Calibrating voltage	$1\text{mV} \pm 2\%$
Voltage tolerance	$\pm 500\text{mV}$
Time constant	$> 3.2\text{s}$
Frequency response	$0.05\text{Hz} \sim 160\text{Hz} (-3\text{db})$
Noise level	$\leq 15 \mu\text{V}_{\text{p-p}}$
Threshold	$\leq 20 \mu\text{V}$
Paper speed	6.25, 12.5, 25, 50mm / st ($\pm 3\%$)
Sensitivity	Auto, 2.5, 5, 10, 20, 40mm/mV
Recorder	High resolution thermal printer
Safety standard	GB9706.1-2007, IEC Class I, Type CF

* Specifications subject to change without notice.



Nidek Medical India Private Limited

Head Office: 238 B, A.J.C. Bose Road, Kolkata 700 020

Tel: (033) 4017 7000, E-mail: sales@nidekindia.com

Website: www.nidekindia.com

Branch Locations:

New Delhi | Mumbai | Chennai | Hyderabad | Lucknow

Pune | Nagpur | Ahmedabad | Chandigarh | Indore

Bangalore | Kochi | Guwahati

Customer Care: +91 93395 66444