

Cardiac Monitoring

NIDEK ECG Machine ECG-703

Portable, Precise, User-friendly

Nidek India's Electrocardiograph machines are one of the basic diagnostic tools for cardiologists and general physicians.

They can be used to perform routine diagnostic checks. The machines are simple to use and precise, thus, making the diagnosis easy to perform and analyze.

Nidek's ECG machines are available in several models suited to the physicians and patient's needs and requirements.

Product descriptions appear in the following pages.

Key Features

- 3 channel ECG machine with interpretation
- Colour display
- Manual / Auto / Rhythm / Exam.
- 80mm 3-ch recording with rhythm lead
- Built in rechargeable battery
- Adapt to 110-230V, 50 / 60Hz AC power supply
- Real-time ECG waveform, cases save / replay
- Built in USB / RS232 interface (ECGNET software optional)
- Pacemaker detection



Specifications

Display	LCD 3.5"
Paper	80mm X 20m roll / fold paper
Dimension	285 X 265 X 72 Unit (mm)

Standard Specifications

Lead	Standard 12 leads
Acquisition	Simultaneous 12 leads
Resolution	12bit / 1000Hz
Operation mode	Manual / Auto / Rhythm / Exam.
Filter	AC, EMG Filter
Drift filter	Anti-Drift System: 0.15 Hz (Adaptive)
CMRR	≥ 120 dB, with AC Filter
Input circuit	Floating, Protection against defibrillator effect
Input circuit current	≥ 0.05 μ A
Input Impedance	≥ 50 M Ω
Patient current leakage	< 10 μ A
Calibrating voltage	1 mV $\pm 5\%$
Sampling rate	4000 Hz / channel
HR Range	30 - 300 BPM (± 1 BPM)
Patient cases storage	400 patient cases. (Internal SD card in optional)
Noise level	≤ 30 μ V _{p-p}
Communication interface	USB, UART
Paper speed	12.5, 25, 50mm / st ($\pm 3\%$)
Sensitivity	2.5, 5, 10, 20, 40mm/mV ($\pm 3\%$)
Recorder	Thermal dot matrix recorder
Safety standard	GB9706.1-2007, IEC Class I, Type CF

* Specifications subject to change without notice.



Nidek Medical India Private Limited

Head Office: 238 B, A.J.C. Bose Road, Kolkata 700 020

Tel: (033) 4017 7000, E-mail: sales@nidekindia.com

Website: www.nidekindia.com

Branch Locations:

New Delhi | Mumbai | Chennai | Hyderabad | Lucknow

Pune | Nagpur | Ahmedabad | Chandigarh | Indore

Bangalore | Kochi | Guwahati

Customer Care: +91 93395 66444