

INTELLIGENT UNIVERSAL VENTILATORS FOR ICU

FLEXIMAG PLUS



Key Features

- Universal Ventilators from extreme-low-weight patients to adults.
- Safe ventilation suggestions as per patient height / weight.
- Touch operated 15 inch display.
- Upgradeable to Capnography or Oximetry with Masimo Sensor (Optional).
- 3 Hour of Battery back up.
- Remote Assistance in Troubleshooting.
- Intelligent Alarm System
- Low Cost of Usage & Ownership.
- Can be used with Nasal Prongs.
- Works with a standard Dual Limb Respiratory Circuit.
- Breath Synchronized Nebulizer
- TGI: Tracheal Gas Insufflation.
- Upgradable to Paramagnetic Oxygen Sensing.

Modes

VCV/VCV-AC, PCV/PCV-AC, PRVC, PLV/PLV with VG, V-SIMV+PS, P-SIMV+PS, DualPAP/APRV, CPAP/PSV, MMV, VS, NIV.

PARAMETERS ADJUSTMENT

Type of patient	Adult, Pediatric and Neonatal
Tidal volume	2 to 3000 ml
Respiratory rate	0 to 200 rpm
Inspiratory flow	1 to 180 L/min
Rise time	0 to 2.0 s
Inspiratory time	0.05 to 30 s
Inspiratory pressure	0 to 120 cmH ₂ O (or hPa or mbar)
Peep	0 to 50 cmH ₂ O (or hPa or mbar)
Support Pressure/ps	0 to 120 cmH ₂ O (or hPa or mbar)
Flow cycling (% of peak flow)	5 to 80 %
Pressure trigger	0.0 to -20 cmH ₂ O (or hPa or mbar)
Flow trigger	0.0 to 30 L/min
Ratio I:E	1:599 to 299:1
O ₂ concentration	21 to 100%
Type of inspiratory flow	Constant, decelerating, accelerating and sine
Inspiratory and expiratory pause	0.1 to 30 s

ALARMS

Minute volume / Total volume	high / low
Respiratory rate	high / low
Inspiratory pressure	high / low
Maximum pressure	high / low
Peep	high / low
FiO ₂	high / low
Apnea time	OFF, 0 to 60 s
Automatic alarm adjustments	OFF, 10%, 20% and 30%

VENTILATION MODES

VCV / VCV-AC; PCV / PCV-AC; PRVC; PLV / PLV with VG; V-SIMV + PS; P-SIMV + PS; DualPAP / APRV; CPAP/PSV; MMV; VS; NIV

OPERATING CONDITIONS SPECIFICATIONS

Electrical power supply	100 to 240 V, 50/60 Hz
12 V _{DC} external	yes (optional)
Battery	210 minutes
O ₂ inlet	29 to 87 psi (200 to 600 kPa)
AR gas inlet	29 to 87 psi (200 to 600 kPa)
Temperature	-10 to 50°C (14 to 122°F)
Barometric pressure	600 to 1100 cmH ₂ O (or hPa or mbar)
Relative humidity	15 to 95%

MONITORING

Curve	PxT, FxT and VxT / SpO ₂ / CO ₂
Loops	PxF, VxF, PxV
Different colors	Insp. and exp. phases, trigger modes and windows
Bargraph	Peak pressure, plateau or instant
FiO ₂	Galvanic or paramagnetic cell (optional)
Optional monitoring	Capnography or Oximetry
Numerical value	Tidal volume and Minute volume; Respiratory rate; Inspiratory and expiratory time; Max and mean airway pressure and plateau pressure; Peep; I:E ratio; FiO ₂ ; Resistance and compliance (static and dynamic); P0.1

USER INTERFACE

Type and Size	TFT-LCD touchscreen 15"
Weight	18,0 kg (39,7 lbs)
Dimensions W x H x D	453 x 1335 x 542mm (17.8 x 52.6 x 21.3 inches)
Communication/Interface	RS-232C ports
Remote Technical Assistance	Remote Assistance Magnamed (ARM)

OTHERS OPERATIONS


Nebulizer	Synchronized with inspiration
Tracheal gas insufflation (TGI)	Synchronized with expiration
Trend	72 hour
Volume compensation - temperature and humidity	Fleximag Plus: BTPS and ATPD
Auxiliary pressure measurement	Using esophageal balloon or pressure at the carina

GENERAL SPECIFICATIONS

Stand by	on/off
Manual cycles	yes
Graphic freeze	yes
Sigh	Fleximag Plus: yes
Flow sensor	Proximal or Distal

MECHANICAL VENTILATION EVALUATION

P0.1	yes
Slow Vital Capacity	yes
PV flex	yes
Plmáx (NIF)	yes
Trapped Volume	yes



MAGNAMED REMOTE ASSISTANCE
TECHNOLOGY MAKES DIFFERENCE

MAGNAMED
Intelligent innovation for life

CE 2460



Nidek Medical India Private Limited
Head Office: 238 B, A.J.C. Bose Road, Kolkata 700 020
Tel: (033) 4017 7000, E-mail: sales@nidekindia.com
Website: www.nidekindia.com

Branch Locations:

New Delhi | Mumbai | Chennai | Hyderabad | Lucknow
Pune | Nagpur | Ahmedabad | Chandigarh | Indore
Bangalore | Kochi | Guwahati

Customer Care: +91 93395 66444