

# INFUSION PUMP

## TERRA

### Features

- Precise and Reliable. Patient safety being of utmost importance. All functions are tested and meet all required guidelines
- Simple - The patient-focussed user interface helps reduce training time by providing therapy related information in an easy to understand format
- Smart - Graphs and reports provide trending data and delivery history for easy patient assessment, continuous quality improvement, and better patient outcomes



## TECHNICAL SPECIFICATIONS

Model	TERRA
Pumping Mechanism	Curvilinear peristaltic
IV Set	Compatible with IV sets of any standard
IV Set Data	8 IV set calibration data storage
Flow Rate	1-1300 ml/h (in 0.1 ml/h increments)
Purge, Bolus	Purge when pump stops, bolus when pump starts, rate at 1100 ml/h
Accuracy	±3%
VTBI	1-9999 ml
Infusion Mode	ml/h, drop/min, time-based
KVO Rate	1-5 ml/h (in 0.1 ml/h increments)
Alarms	Occlusion, air-in-line, door open, end program, low battery, end battery, AC power off, motor malfunction, system malfunction, reminder alarm
Additional Features switching,	Real-time infused volume, automatic power mute key, purge, bolus, system memory, key locker, nurse call
Occlusion Sensitivity	5 levels
Air-in-line Detection	Ultrasonic detector
Nurse Call	Available
Wireless management	Optional
Drop Sensor	Optional
Power Supply, AC	110/230 V (optional), 50-60 Hz, 20 VA
Battery	9.6±1.6 V, rechargeable
Battery Life	6 hours at 30 ml/h
Size	233*146*269 mm
Weight	3 kg
Safety Classification	Class I, type CF, CE



Nidek Medical India Private Limited

Head Office: 238 B, A.J.C. Bose Road, Kolkata 700 020

Tel: (033) 4017 7000, E-mail: sales@nidekindia.com

Website: www.nidekindia.com

Branch Locations:

New Delhi | Mumbai | Chennai | Hyderabad | Lucknow

Pune | Nagpur | Ahmedabad | Chandigarh | Indore

Bangalore | Kochi | Guwahati

Customer Care: +91 93395 66444

