




FINGERTIP PULSE OXIMETER

MODEL-6900


Hand Held
Plethysmography with Graph
Pulse Rate

Assessing a patient's need for oxygen is the most essential element to life; no human life thrives in the absence of oxygen. Nidek Fingertip oximeter is useful in any setting where a patient's oxygenation is unstable, including intensive care, operating, recovery, emergency and hospital ward settings, or at home, for assessment of any patient's oxygenation



SPECIFICATIONS

SPO2, PR, SPO2 bar graph and pulse waveform display
Accuracy: SPO2 \pm 2%; Pulse \pm 2 BPM
Range: SPO2 - 70-99%; Pulse 30-235 BPM
Dual color OLED display with 4 direction view
PI index to show pulsating blood situation
Power supply: 2 pcs standard AAA batteries
Batterly level indication
Better shading level. Hingh accuracy and durability
Operation menu for function setting
Automatic switch off function for power saving
Suitable to Neonatal, Paediatric and Adult use
Weight: 50 grams (including two AAA batteries)



Nidek Medical India Private Limited

Head Office: 238 B, A.J.C. Bose Road, Kolkata 700 020

Tel: (033) 4017 7000, E-mail: sales@nidekindia.com

Website: www.nidekindia.com

Branch Locations:

New Delhi | Mumbai | Chennai | Hyderabad | Lucknow

Pune | Nagpur | Ahmedabad | Chandigarh | Indore

Customer Care: **+91 93395 66444**